


Product Information	
Catalog No.	RA12409
Product name	Anti-Human Lp-PLA2/PLA2G7/PAF-AH Antibody (SAA1765)
Applications	ELISA
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	SAA1765
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Platelet-activating factor acetylhydrolase, PAFAH, PAF 2-acylhydrolase, PAF acetylhydrolase, 2-acetyl-1-alkylglycerophosphocholine esterase, 1-alkyl-2-acetyl-glycerophosphocholine esterase, LDL-PLA(2), LDL-associated phospholipase A2, PLA2G7, Group-VIIA phospholipase A2, gVIIA-PLA2
Accession	Q13093
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human Lp-PLA2/PLA2G7/PAF-AH Antibody (SAA1765)

SDS-PAGE