


Product Information	
Catalog No.	RA12456
Product name	Anti-Human PTGER4/EP4 Receptor Antibody (SAA1812)
Applications	ELISA
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	SAA1812
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	PTGER2, Prostaglandin E2 receptor EP4 subtype, PTGER4, PGE receptor EP4 subtype, Prostanoid EP4 receptor, PGE2 receptor EP4 subtype
Accession	P35408
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human PTGER4/EP4 Receptor Antibody (SAA1812)

SDS-PAGE