


Product Information	
Catalog No.	RA12483
Product name	Anti-Human ST14 Antibody (SAA1839)
Applications	ELISA
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	SAA1839
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Tumor-associated differentially-expressed gene 15 protein, Matriptase, Membrane-type serine protease 1, SNC19, TADG15, Serine protease 14, MT-SP1, Prostamin, ST14, Serine protease TADG-15, Suppressor of tumorigenicity 14 protein, PRSS14
Accession	Q9Y5Y6
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human ST14 Antibody (SAA1839)

SDS-PAGE