


Product Information	
Catalog No.	RA12526
Product name	Anti-Human NGLY1 Antibody (SAA1882)
Applications	ELISA
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	SAA1882
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	PNGase, N-glycanase 1, NGLY1, Peptide-N(4)-(N-acetyl-beta-glucosaminyl)asparagine amidase, hPNGase, PNG1, Peptide:N-glycanase
Accession	Q96IV0
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human NGLY1 Antibody (SAA1882)

SDS-PAGE