


Product Information	
Catalog No.	RA12556
Product name	Anti-Human TMEM106B Antibody (SAA1912)
Applications	ELISA, FCM, SPR
Species reactivity	Human, Mouse
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	SAA1912
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	TMEM106B, Transmembrane protein 106B
Accession	Q9NUM4, Q80X71
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human TMEM106B Antibody (SAA1912)

SDS-PAGE