


Product Information	
Catalog No.	RA12618
Product name	Anti-Human SLC4A7 Antibody (SAA1974)
Applications	ELISA
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	VHH-mFc
Clone ID	SAA1974
Form	Liquid
Storage buffer	30 mM Histidine, pH5.8, 10% sucrose, 0.02% Tween80
Concentration	0.44 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Electroneutral Na/HCO ₃ cotransporter, Sodium bicarbonate cotransporter 3, NBC3, Sodium bicarbonate cotransporter 2, NBCn1, Solute carrier family 4 member 7, Bicarbonate transporter, SLC4A6, NBC2, SBC2, Sodium bicarbonate cotransporter 2b, NBC2B, BT, SLC4A7
Accession	Q9Y6M7
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human SLC4A7 Antibody (SAA1974)

SDS-PAGE