


Product Information	
Catalog No.	RA12672
Product name	Anti-Neisseria meningitidis NadA/Neisseria adhesin A Antibody (G9)
Applications	DB, ELISA, WB
Species reactivity	Neisseria meningitidis serogroup B (strain MC58)
Host species	Human
Clonality	Monoclonal
Isotype	VHH-hFc
Clone ID	G9
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	nadA, Neisseria adhesin A
Accession	Q9JXK7
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000; WB: 1:500-1:2000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Neisseria meningitidis NadA/Neisseria adhesin A Antibody (G9)

SDS-PAGE