


Product Information	
Catalog No.	RA12756
Product name	Anti-Human HLA-G Antibody (38426#)
Applications	FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG1
Clone ID	38426#
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	HLA-G, HLA-G, sHLA-G, HLA G antigen, MHC class I antigen G, HLA-6.0, HLA class I histocompatibility antigen, alpha chain G
Accession	P17693
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human HLA-G Antibody (38426#)

SDS-PAGE