


Product Information


Catalog No.	RA12802
Product name	Anti-Human IL13 Antibody (SAA2016)
Applications	ELISA, WB
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG4, kappa
Clone ID	SAA2016
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	IL13, NC30, Interleukin-13, IL-13
Accession	P35225
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000; WB: 1:500-1:2000
Note	For research use only.

Image Data




SEC-HPLC detection for Anti-Human IL13 Antibody (ABS1981).

Bioactivity



SDS-PAGE for Anti-Human IL13 Antibody (SAA2016).

SDS-PAGE



Western blot analysis was performed using anti-IL13 monoclonal antibody at 1ug/ml on various samples.

Lane 1: recombinant human IL13 (Catalog No: HW632012), 50ng

Lane 2: recombinant human IL13 (Catalog No: HW632012), 30ng

Lane 3: recombinant human IL13 (Catalog No: HW632012), 20ng

Western blot