


Product Information	
Catalog No.	RA12849
Product name	Anti-SAPAP1 Antibody (N238/30)
Applications	IHC, WB
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	N238/30
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	DAP-1, DAP1, SAP90/PSD-95-associated protein 1, PSD-95/SAP90-binding protein 1, Guanylate kinase-associated protein, hGKAP, GKAP, Disks large-associated protein 1, SAPAP1, DLGAP1
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	WB: 1:500-1:2000
Note	For research use only.

Image Data



SDS-PAGE for Anti-SAPAP1 Antibody (N238/30)

SDS-PAGE