


Product Information	
Catalog No.	RA12873
Product name	Anti-WRP/SGAP3 Antibody (N285/28)
Applications	IHC, WB
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	N285/28
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	SGAP3, srGAP3, SLIT-ROBO Rho GTPase-activating protein 3, WRP, WAVE-associated Rac GTPase-activating protein, KIAA0411, Rho GTPase-activating protein 14, SGAP2, Mental disorder-associated GAP, ARHGAP14, KIAA1156, MEGAP
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	WB: 1:500-1:2000
Note	For research use only.

Image Data



SDS-PAGE for Anti-WRP/SGAP3 Antibody (N285/28)

SDS-PAGE