


Product Information	
Catalog No.	RA12979
Product name	Anti-Gs protein, alpha subunit Antibody (N192/24)
Applications	IHC, WB
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	N192/24
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	POLA1, POLA, DNA polymerase alpha catalytic subunit p180, DNA polymerase alpha catalytic subunit
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	WB: 1:500-1:2000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Gs protein, alpha subunit Antibody (N192/24)

SDS-PAGE