


Product Information	
Catalog No.	RA13010
Product name	Anti-Kv1.2 K ⁺ channel Antibody (L76/44)
Applications	IHC, WB
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	L76/44
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Voltage-gated K(+) channel HuKIV, Voltage-gated potassium channel subunit Kv1.2, NGK1, KCNA2, Potassium voltage-gated channel subfamily A member 2, Voltage-gated potassium channel HBK5
Accession	P16389
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	WB: 1:500-1:2000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Kv1.2 K⁺ channel Antibody (L76/44)

SDS-PAGE