


| Product Information       |                                                                                                                                                              |
|---------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Catalog No.               | RA13100                                                                                                                                                      |
| Product name              | Anti-GalR1 Antibody (N371B/37)                                                                                                                               |
| Applications              | IHC, WB                                                                                                                                                      |
| Host species              | Mouse                                                                                                                                                        |
| Clonality                 | Monoclonal                                                                                                                                                   |
| Isotype                   | IgG2a, kappa                                                                                                                                                 |
| Clone ID                  | N371B/37                                                                                                                                                     |
| Form                      | Liquid                                                                                                                                                       |
| Storage buffer            | 0.01M PBS, pH 7.4.                                                                                                                                           |
| Concentration             | 1 mg/ml                                                                                                                                                      |
| Purity                    | >95% as determined by SDS-PAGE.                                                                                                                              |
| Purification              | Protein A/G purified from cell culture supernatant.                                                                                                          |
| Target                    | GAL1-R, Galanin receptor type 1, GALR-1, GALR1, GALNR1, GALNR                                                                                                |
| Endotoxin level           | Please contact with the lab for this information.                                                                                                            |
| Stability and Storage     | Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term. |
| Product Usage Information | WB: 1:500-1:2000                                                                                                                                             |
| Note                      | For research use only.                                                                                                                                       |

## Image Data



SDS-PAGE for Anti-GalR1 Antibody (N371B/37)

SDS-PAGE