


Product Information	
Catalog No.	RA13226
Product name	Anti-GluN2B/NR2B glutamate receptor Antibody (N59/36)
Applications	IHC, WB
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	N59/36
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	GluN2B, Glutamate receptor ionotropic, NMDA 2B, GRIN2B, NMDAR2B, NR2B, Glutamate [NMDA] receptor subunit epsilon-2, N-methyl-D-aspartate receptor subunit 3, N-methyl D-aspartate receptor subtype 2B, hNR3, NR3
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	WB: 1:500-1:2000
Note	For research use only.

Image Data



SDS-PAGE for Anti-GluN2B/NR2B glutamate receptor Antibody (N59/36)

SDS-PAGE