


Product Information	
Catalog No.	RA13233
Product name	Anti-JIP-2/IB-2 Antibody (N135/13)
Applications	IHC, WB
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	N135/13
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	MAPK8IP2, Islet-brain-2, IB-2, JNK MAP kinase scaffold protein 2, C-Jun-amino-terminal kinase-interacting protein 2, JIP2, JNK-interacting protein 2, Mitogen-activated protein kinase 8-interacting protein 2, PRKM8IPL, IB2, JIP-2
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	WB: 1:500-1:2000
Note	For research use only.

Image Data



SDS-PAGE for Anti-JIP-2/IB-2 Antibody (N135/13)

SDS-PAGE