


Product Information	
Catalog No.	RA13302
Product name	Anti-Cav2.3 Ca ²⁺ channel Antibody (N493/92)
Applications	IHC, WB
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	N493/92
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	CACH6, CACNA1E, BII, Calcium channel, L type, alpha-1 polypeptide, isoform 6, Voltage-dependent R-type calcium channel subunit alpha-1E, Brain calcium channel II, CACNL1A6, Voltage-gated calcium channel subunit alpha Cav2.3
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	WB: 1:500-1:2000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Cav2.3 Ca²⁺ channel Antibody (N493/92)

SDS-PAGE