


Product Information	
Catalog No.	RA13312
Product name	Anti-Frataxin Antibody (N191/7)
Applications	IHC, WB
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	N191/7
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	FXN, Frataxin(81-210), Friedreich ataxia protein, m56-FXN, Frataxin, mitochondrial, Fxn, d-FXN, m78-FXN, i-FXN, m81-FXN, X25, FRDA
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	WB: 1:500-1:2000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Frataxin Antibody (N191/7)

SDS-PAGE