


Product Information	
Catalog No.	RA13358
Product name	Anti-Histone H3.3 Antibody (N158/28)
Applications	IHC, WB, ELISA
Species reactivity	Human, Mouse, Rat
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	N158/28
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Histone H3.3C, H3-5, Histone H3.5, H3F3C
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000; WB: 1:500-1:2000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Histone H3.3 Antibody (N158/28)

SDS-PAGE