


Product Information	
Catalog No.	RA13410
Product name	Anti-TRIP8b Antibody (N212/17)
Applications	IHC, WB
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	N212/17
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Tetratricopeptide repeat-containing Rab8b-interacting protein, PEX5-related protein, PEX5-like protein, Peroxisome biogenesis factor 5-like, Pex5Rp, Peroxin-5-related protein, TRIP8b, PXR2, PEX5R, PEX2-related protein, PEX5L
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	WB: 1:500-1:2000
Note	For research use only.

Image Data



SDS-PAGE for Anti-TRIP8b Antibody (N212/17)

SDS-PAGE