


Product Information	
Catalog No.	RA13555
Product name	Anti-P2X purinoceptor 7 Antibody (N256/43)
Applications	IHC, WB
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	N256/43
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	P2RX7, Purinergic receptor, P2X purinoceptor 7, P2X7, P2Z receptor, ATP receptor
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	WB: 1:500-1:2000
Note	For research use only.

Image Data



SDS-PAGE for Anti-P2X purinoceptor 7 Antibody (N256/43)

SDS-PAGE