


Product Information	
Catalog No.	RA13605
Product name	Anti-SNAT1 Antibody (N104/18)
Applications	IHC, WB
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	N104/18
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Sodium-coupled neutral amino acid transporter 1, Solute carrier family 38 member 1, SLC38A1, N-system amino acid transporter 2, SNAT1, ATA1, NAT2, System A amino acid transporter 1, SAT1, System N amino acid transporter 1, Amino acid transporter A1
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	WB: 1:500-1:2000
Note	For research use only.

Image Data



SDS-PAGE for Anti-SNAT1 Antibody (N104/18)

SDS-PAGE