


Product Information	
Catalog No.	RA13687
Product name	Anti-KChIP1 K+ channel Antibody (K59/19)
Applications	IHC, WB
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	K59/19
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	A-type potassium channel modulatory protein 1, Vesicle APC-binding protein, KCNIP1, Potassium channel-interacting protein 1, KChIP1, Kv channel-interacting protein 1, VABP, KCHIP1
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	WB: 1:500-1:2000
Note	For research use only.

Image Data



SDS-PAGE for Anti-KChIP1 K+ channel Antibody (K59/19)

SDS-PAGE