


Product Information	
Catalog No.	RA13706
Product name	Anti-DNAH1 Antibody (N251/14)
Applications	IHC, WB
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	N251/14
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	hDHC7, Heat shock regulated protein 1, HSRF-1, Axonemal beta dynein heavy chain 1, DHC7, KIAA1410, DNAH1, Ciliary dynein heavy chain 1, DNAHC1, Dynein heavy chain 1, axonemal
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	WB: 1:500-1:2000
Note	For research use only.

Image Data



SDS-PAGE for Anti-DNAH1 Antibody (N251/14)

SDS-PAGE