


Product Information	
Catalog No.	RA13770
Product name	Anti-IRSp53/BAIAP2 Antibody (L117/1)
Applications	IHC, WB
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	L117/1
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	BAI-associated protein 2, Insulin receptor substrate p53, IRS-58, BAI1-associated protein 2, IRSp53, Insulin receptor substrate p53/p58, Fas ligand-associated factor 3, BAIAP2, FLAF3, Insulin receptor substrate protein of 53 kDa, Brain-specific angiogenesis inhibitor 1-associated protein 2, Protein BAP2, IRSp53/58
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	WB: 1:500-1:2000
Note	For research use only.

Image Data



SDS-PAGE for Anti-IRSp53/BAIAP2 Antibody (L117/1)

SDS-PAGE