


Product Information	
Catalog No.	RA13791
Product name	Anti-MMACHC Antibody (N230/8)
Applications	IHC, WB
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	N230/8
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	MMACHC, Methylmalonic aciduria and homocystinuria type C protein, Cyanocobalamin reductase (cyanide-eliminating), Alkylcobalamin:glutathione S-alkyltransferase, Cyanocobalamin reductase / alkylcobalamin dealkylase, CblC
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	WB: 1:500-1:2000
Note	For research use only.

Image Data



SDS-PAGE for Anti-MMACHC Antibody (N230/8)

SDS-PAGE