


Product Information	
Catalog No.	RA13823
Product name	Anti-Human MAP3K2 Antibody (E6#)
Applications	ELISA, SPR, WB
Species reactivity	Homo sapiens (Human), Mus musculus (Mouse)
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Clone ID	E6#
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	MAP3K2, MAPKKK2, MEKK2, Mitogen-activated protein kinase kinase kinase 2, MAPK/ERK kinase kinase 2, MEK kinase 2, MEKK 2
Accession	Q9Y2U5, Q61083
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000; WB: 1:500-1:2000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human MAP3K2 Antibody (E6#)

SDS-PAGE