


Product Information	
Catalog No.	RA13834
Product name	Anti-Optineurin Antibody (6A7.1)
Applications	IHC, WB
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Clone ID	6A7.1
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	NRP, TFIIIA-IntP, Huntingtin-interacting protein 7, HIP-7, Huntingtin yeast partner L, Transcription factor IIIA-interacting protein, NEMO-related protein, E3-14.7K-interacting protein, GLC1E, HIP7, HYPL, FIP-2, Huntingtin-interacting protein L, Optic neuropathy-inducing protein, OPTN, Optineurin, FIP2
Accession	Q96CV9
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	WB: 1:500-1:2000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Optineurin Antibody (6A7.1)

SDS-PAGE