


Product Information	
Catalog No.	RA13886
Product name	Anti-Bacteriophage 933W stxB2/SLT-2b/VT2 Nanobody (Nb113)
Applications	ELISA, SPR, WB
Species reactivity	Escherichia phage 933W (Bacteriophage 933W)
Host species	Alpaca
Clonality	Monoclonal
Isotype	VHH-8His-Cys-tag
Clone ID	Nb113
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Purified by Nickel column.
Target	Shiga-like toxin 2 subunit B, SLT-2 B subunit, SLT-2b, SLT-IIb, Verocytotoxin 2 subunit B, Verotoxin 2 subunit B, stxB2, stx2B, L0104
Accession	P09386
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000; WB: 1:500-1:2000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Bacteriophage 933W stxB2/SLT-2b/VT2 Nanobody (Nb113)

SDS-PAGE