


Product Information	
Catalog No.	RA13951
Product name	InVivoMAb Anti-Human APOE4 Antibody (HJ15.1)
Applications	ELISA, FuncS
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	HJ15.1
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	0.39 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Apo-E, Apolipoprotein E, APOE
Accession	P02649
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for InVivoMAb Anti-Human APOE4 Antibody (HJ15.1)

SDS-PAGE