

Product Information	
Catalog No.	RA13956
Product name	Anti-SARS-CoV-2 Spike protein IgM Antibody (S309)
Applications	ELISA, Neutralization
Species reactivity	SARS-CoV-2(2019-nCoV)
Host species	Human
Clonality	Monoclonal
Isotype	IgM
Clone ID	S309
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein L purified from cell culture supernatant.
Target	Receptor-Binding Domain, RBD, Spike glycoprotein, S glycoprotein, E2, Peplomer protein, Spike protein S1, S, SARS-CoV-2, 2019-nCoV
Accession	P0DTC2
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.