


Product Information	
Catalog No.	RA13965
Product name	Anti-Human ROR1/NTRKR1 Antibody (4A5)
Applications	ELISA, FCM
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Clone ID	4A5
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Neurotrophic tyrosine kinase, receptor-related 1, Inactive tyrosine-protein kinase transmembrane receptor ROR1, NTRKR1, ROR1
Accession	Q01973
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human ROR1/NTRKR1 Antibody (4A5)

SDS-PAGE