


Product Information	
Catalog No.	RA13992
Product name	Anti-DENV-2 Envelope protein E/EDIII A strand Antibody (E76#)
Applications	ELISA
Species reactivity	Dengue virus type 2 (strain Puerto Rico/PR159-S1/1969) (DENV-2)
Host species	Human
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	E76#
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Genome polyprotein, Envelope protein E
Accession	P12823
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-DENV-2 Envelope protein E/EDIII A strand Antibody (E76#)

SDS-PAGE