


Product Information	
Catalog No.	RA14038
Product name	InVivoMAb Anti-Human HLA-A*11:01 Antibody (2E3#)
Applications	ELISA, FCM, FuncS
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	2E3#
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	HLAA, Human leukocyte antigen A, HLA class I histocompatibility antigen, A alpha chain, HLA-A, HLA-A*11:01
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for InVivoMAb Anti-Human HLA-A*11:01 Antibody (2E3#).

SDS-PAGE