


Product Information	
Catalog No.	RA14132
Product name	Anti-Human ZSCAN23 Antibody (C482)
Applications	ELISA, IF, IP
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	C482
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	ZSCAN23, ZNF390, ZNF453, Zinc finger and SCAN domain-containing protein 23, Zinc finger protein 390, Zinc finger protein 453
Accession	Q3MJ62
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human ZSCAN23 Antibody (C482)

SDS-PAGE