


Product Information	
Catalog No.	RA14151
Product name	Anti-Human ZBTB21 Antibody (C85)
Applications	ELISA, IF, IP
Species reactivity	Human
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	C85
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	ZBTB21, KIAA1227, ZNF295, Zinc finger and BTB domain-containing protein 21, Zinc finger protein 295
Accession	Q9ULJ3
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human ZBTB21 Antibody (C85)

SDS-PAGE