


Product Information	
Catalog No.	RA14169
Product name	InVivoMAb Anti-SARS-CoV-2 RBD Antibody (PR961)
Applications	ELISA, Neutralization
Species reactivity	SARS-CoV-2
Host species	Human
Clonality	Monoclonal
Isotype	IgG1
Clone ID	PR961
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Spike glycoprotein, S glycoprotein, E2, Peplomer protein, RBD, Spike protein S1, S, Receptor-Binding Domain, S 6P trimer, SARS-CoV-2
Accession	P0DTC2
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for InVivoMAb Anti-SARS-CoV-2 RBD Antibody (PR961)

SDS-PAGE