

Product Information	
Catalog No.	RA14192
Product name	Anti-Human CD16b/FCGR3B Nanobody (C21)
Applications	ELISA, WB, IF, FuncS
Species reactivity	Human
Host species	Alpaca
Clonality	Monoclonal
Isotype	VHH-8His-Cys-tag
Clone ID	C21
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Purified by Nickel column.
Target	IGFR3, FCGR3B, CD16b, Low affinity immunoglobulin gamma Fc region receptor III-B, Fc-gamma RIII, IgG Fc receptor III-1, FCGR3, CD16B, Fc-gamma RIII-beta, Fc-gamma RIIIb, FCG3, FcRIII, FcR-10, FcRIIIb
Accession	O75015
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000; WB: 1:500-1:2000
Note	For research use only.