

Product Information	
Catalog No.	RA14262
Product name	Anti-Bovine Serum Albumin/BSA Antibody (12D6), HRP
Applications	ELISA (Det), WB
Species reactivity	Bovine
Host species	Human
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	12D6
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4, 50% Glycerol, 0.05% Proclin 300.
Concentration	0.39 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Albumin, BSA, Bos d 6, ALB
Accession	P02769
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000; WB: 1:500-1:2000
Note	For research use only.