


Product Information	
Catalog No.	RA14296
Product name	InVivoMAb Anti-Mouse PD-L1 & CD3e Bispecific Antibody (Iv0242)
Applications	FuncS
Species reactivity	Mouse
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG
Clone ID	Iv0242
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	0.8 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Bispecific, Programmed cell death 1 ligand 1, PD-L1, PDCD1 ligand 1, Programmed death ligand 1, B7 homolog 1, B7-H1, CD274, Cd274, B7h1, Pdcd1l1, Pdcd1lg1, Pdl1, T-cell surface glycoprotein CD3 epsilon chain, T-cell surface antigen T3/Leu-4 epsilon chain, CD3e, Cd3e
Accession	Q9EP73 & P22646
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for InVivoMAb Anti-Mouse PD-L1 & CD3e Bispecific Antibody (lv0242)

SDS-PAGE



Western blot analysis was performed using anti-/ monoclonal antibody at 1ug/mL on various samples.

Lane 1: ATGE3052 200ng

Lane 2: 1* loading buffer

Western blot