


Product Information	
Catalog No.	RA14299
Product name	InVivoMAb Anti-Mouse PD-L1 & CD28 Bispecific Antibody (Iv0245)
Applications	FuncS
Species reactivity	Mouse
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG
Clone ID	Iv0245
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Bispecific, Programmed cell death 1 ligand 1, PD-L1, PDCD1 ligand 1, Programmed death ligand 1, B7 homolog 1, B7-H1, CD274, Cd274, B7h1, Pdcd1l1, Pdcd1lg1, Pdl1, T-cell-specific surface glycoprotein CD28, CD28, Cd28
Accession	Q9EP73 & P31041
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for InVivoMAb Anti-Mouse PD-L1 & CD28 Bispecific Antibody (lv0245)

SDS-PAGE