


Product Information	
Catalog No.	RA14314
Product name	Anti-Helicobacter pylori vacA Antibody (N03-F10)
Applications	ELISA, FCM, SPR, WB
Species reactivity	Helicobacter pylori (Campylobacter pylori), Helicobacter pylori (strain ATCC 700392 / 26695) (Campylobacter pylori)
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	N03-F10
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Vacuolating cytotoxin autotransporter, Vacuolating cytotoxin, Vacuolating cytotoxin translocator, vacA, Vacuolating cytotoxin autotransporter, Vacuolating cytotoxin, Vacuolating cytotoxin translocator, vacA, HP_0887
Accession	Q48245, P55981
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000; WB: 1:500-1:2000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Helicobacter pylori vacA Antibody (N03-F10)

SDS-PAGE