


Product Information	
Catalog No.	RA14331
Product name	Anti-Human CD112/Nectin-2/PVRL2 Nanobody (SAA2117)
Applications	ELISA
Species reactivity	Human
Host species	Alpaca
Clonality	Monoclonal
Isotype	VHH-8His-Cys-tag
Clone ID	SAA2117
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Purified by Nickel column.
Target	NECTIN2, HVEB, PRR2, PVRL2, Nectin-2, CD112, Herpes virus entry mediator B, Herpesvirus entry mediator B, HveB, Nectin cell adhesion molecule 2 receptor-related protein 2
Accession	Q92692
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.

Image Data



SDS-PAGE for Anti-Human CD112/Nectin-2/PVRL2 Nanobody (SAA2117).

SDS-PAGE