

Product Information	
Catalog No.	RA14410
Product name	Anti-Human CD85k/LILRB4 Antibody (h128-3)
Applications	ELISA, FCM
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	h128-3
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Leukocyte immunoglobulin-like receptor subfamily B member 4, B4, CD85 antigen-like family member K, Immunoglobulin-like transcript 3, ILT-3, Leukocyte immunoglobulin-like receptor 5, LIR-5, Monocyte inhibitory receptor HM18, CD85k, LILRB4, ILT3, LIR5
Accession	Q8NHJ6
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.