


Product Information	
Catalog No.	RA14419
Product name	InVivoMAb Anti-Mouse MSLN/Mesothelin Antibody (Ab237)
Applications	ELISA, FCM, Neutralization
Species reactivity	Mouse
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	Ab237
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	Mesothelin, Pre-pro-megakaryocyte-potentiating factor, Megakaryocyte-potentiating factor, MPF, Mesothelin, cleaved form, Msln, Mes, Mpf
Accession	Q61468
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for InVivoMAb Anti-Mouse MSLN/Mesothelin Antibody (Ab237)

SDS-PAGE