

Product Information	
Catalog No.	RA14492
Product name	Anti-Ganglioside GM1 hIgM Antibody (WO1)
Applications	ELISA
Species reactivity	Human
Host species	Human
Clonality	Monoclonal
Isotype	IgM
Clone ID	WO1
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	0.74 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein L purified from cell culture supernatant.
Target	Ganglioside GM1, Ganglioside A2, Ganglioside G4, Ganglioside GGtet1, Ganglioside GI, Monosialoganglioside GM1
Accession	ChEBI: 61048
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Product Usage Information	ELISA: 1:2000-1:20000
Note	For research use only.