


Product Information

Catalog No.	RA14494
Product name	InVivoMAb Anti-HRSV/hMPV Fusion glycoprotein Neutralizing Antibody (M2B6)
Applications	ELISA, Neutralization
Species reactivity	Human respiratory syncytial virus (HRSV), Human metapneumovirus (HMPV)
Host species	Human
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	M2B6
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	F, Fusion glycoprotein F
Accession	P03420, G3KCK8
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only. Not suitable for clinical or therapeutic use.


Image Data

Anti-HRSV/hMPV Fusion glycoprotein Neutralizing Antibody (M2B6) binds with HRSV Fusion glycoprotein F0



Bioactivity

Detects Recombinant HRSV-A2 Pre-F/Fusion glycoprotein F0 Protein, C-His-Strep (Cat No.: VK421011) in indirect ELISA.



Western blot

Western blot analysis was performed using anti-fusion glycoprotein monoclonal antibody at 1ug/mL on various samples.
Lane 1: recombinant human metapneumovirus fusion glycoprotein (Catalog No: VK430011)
Lane 2: negative control