


Product Information	
Catalog No.	RA14497
Product name	InVivoMAb Anti-Mouse CD268/TNFRSF13C Antibody (Iv0266)
Applications	Antagonist, ELISA, Neutralization
Species reactivity	Mouse
Host species	Mouse
Clonality	Monoclonal
Isotype	IgG2a, kappa
Clone ID	Iv0266
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	BR3, TNFRSF13C, BAFF receptor, BLyS receptor 3, Tumor necrosis factor receptor superfamily member 13C, CD268, B-cell-activating factor receptor, BAFFR, BAFF-R
Accession	Q9D8D0
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only. Not suitable for clinical or therapeutic use.

Image Data



SDS-PAGE for InVivoMAb Anti-Mouse CD268/TNFRSF13C Antibody (lv0266)

SDS-PAGE