

Product Information	
Catalog No.	RA14516
Product name	InVivoMAb Anti-Human EDB & CD3E Bispecific Antibody (Iv0268)
Applications	Neutralization
Species reactivity	Human
Clonality	Monoclonal
Isotype	BiTE-bispecific antibody
Clone ID	Iv0268
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Purified by Nickel column.
Target	Bispecific, EDB, EDB-FN, Fibronectin, Cold-insoluble globulin, CIG, T3E, T-cell surface antigen T3/Leu-4 epsilon chain, CD3e, CD3E, T-cell surface glycoprotein CD3 epsilon chain
Accession	P02751 & P07766
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only. Not suitable for clinical or therapeutic use.