


Product Information	
Catalog No.	RA14527
Product name	InVivoMAb Anti-Dog TNFa/TNF-alpha Antibody (Iv0271)
Applications	Blocking, ELISA, FCM
Species reactivity	Dog
Host species	Dog
Clonality	Monoclonal
Isotype	IgG1, kappa
Clone ID	Iv0271
Form	Liquid
Storage buffer	0.01M PBS, pH 7.4.
Concentration	1 mg/ml
Purity	>95% as determined by SDS-PAGE.
Purification	Protein A/G purified from cell culture supernatant.
Target	TNF, Tumor necrosis factor ligand superfamily member 2, N-terminal fragment, ICD2, NTF, TNF-a, TNF-alpha, Tumor necrosis factor, TNFSF2, TNFA, Cachectin, ICD1
Accession	P51742
Endotoxin level	Please contact with the lab for this information.
Stability and Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. Store at 4°C short term (1-2 weeks). Store at -20°C 12 months. Store at -80°C long term.
Note	For research use only. Not suitable for clinical or therapeutic use.


Image Data



SDS-PAGE for InVivoMAb Anti-Dog TNFa/TNF-alpha Antibody (lv0271)

SDS-PAGE

Anti-Dog TNFa/TNF-alpha Antibody (lv0271)
binds with Dog TNFa/TNF-alpha



Detects recombinant dog TNFa/TNF-alpha (Catalog No: CF879011) in indirect ELISA